


Combining Kinetic Ligand Binding and 3D Tumor Invasion Technologies to Assess Drug Residence Time and Anti-metastatic Effects of CXCR4 Inhibitors



Brad Larson and Leonie Rieger, BioTek Instruments, Inc., Winooski, VT
Nicolas Pierre, Cisbio US, Inc., Bedford, MA
Hilary Sherman, Corning Incorporated, Life Sciences, Kennebunk, ME



Introduction

Metastasis, the spread of cancer cells from the original tumor to secondary locations within the body, is linked to approximately 90% of cancer deaths¹. The expression of chemokine receptors, such as CXCR4 and CCR7, is tightly correlated with the metastatic properties of breast cancer cells. *In vivo*, neutralizing the interaction of CXCR4 and its known ligand, SDF1- α (CXCL12), significantly impaired the metastasis of breast cancer cells and cell migration². Traditionally, the discovery of novel agents has been guided by the affinity of the ligand for the receptor under equilibrium conditions, largely ignoring the kinetic aspects of the ligand-receptor interaction. However, awareness of the importance of binding kinetics has started to increase due to accumulating evidence^{3, 4, 5, 6} suggesting that the *in vivo* effectiveness of ligands may be attributed to the time a particular ligand binds to its receptor (drug-target residence time).

Similarly, appropriate *in vitro* cell models have also been lacking to accurately assess the ability of novel therapies to inhibit tumor invasion. Tumors *in vivo* exist as a three-dimensional (3D) mass of multiple cell types, including cancer and stromal cells⁷. Therefore, incorporating a 3D spheroid-type cellular structure that includes co-cultured cell types forming a tumoroid, provides a more predictive model than the use of individual cancer cells cultured on the bottom of a well in traditional two-dimensional (2D) format.

Here we examine the drug-target residence time of various CXCR4 inhibitors using a direct, homogeneous ligand binding assay and CXCR4 expressing cell line in a kinetic format. This inhibitor panel was further tested in a 3D tumor invasion assay to determine whether there is a correlation between the molecule's CXCR4 residence time and inhibition of the phenotypic effect of tumor invasion. MDA-MB-231 breast adenocarcinoma cells, known to be invasive, and metastasize to lung from primary mammary fat pad tumors⁸, were included, in addition to primary human dermal fibroblasts. Cellular analysis algorithms provided accurate quantification

of changes to the original tumoroid structure, as well as invadopodia development. The combination presents an accurate, yet easy-to-use method to assess target-based and phenotypic effects of new, potential anti-metastatic drugs.

Materials and Methods

Materials

Cells

Cryopreserved Tag-lite[®] HEK293 cells stably expressing CXCR4 were provided by Cisbio Assays (Bedford, MA). MDA-MB-231 cells expressing GFP (Catalog No. AKR-211) were purchased from Cell Biolabs, Inc. (San Diego, CA). Human neonatal dermal fibroblasts expressing RFP (Catalog No. cAP-0008RFP) were purchased from Angio-Proteomie (Boston, MA). Both cell types were propagated in Advanced DMEM Medium (Catalog No. 12491-015) plus Fetal Bovine Serum (FBS), 10% (Catalog No. 10437-028) and Pen-Strep-Glutamine, 1x (Catalog No. 10378-016) each from Life Technologies (Carlsbad, CA).

Experimental Components

SDF1- α -d2 and all labeling reagents used for characterization of inhibitor kinetic binding properties were provided by Cisbio Assays. Recombinant Human CXCL12/SDF-1 α (Catalog No. 350-NS-010) was purchased from R&D Systems (Minneapolis, MN). DMEM/F12, HEPES, No Phenol Red Medium (Catalog No. 11039-021) was purchased from Life Technologies. Matrigel Basement Membrane Matrix, Phenol Red-Free (Catalog No. 356237) was purchased from Corning Life Sciences (Corning, NY). AMD 3100 (Catalog No. 3299), AMD 3465 (Catalog No. 4179), WZ 811 (Catalog No. 3951), FC 131 (Catalog No. 4320), IT1t (Catalog No. 4596), and CTCE 9908 (Catalog No. 5130) were purchased from R&D Systems (Minneapolis, MN). 384-well low volume, white, round bottom, non-

Key Words:

3D
3D Cell Culture
Spheroid
Tumoroid
ULA
Drug-Target Residence Time
Residence Time
Ligand Binding
GPCR
CXCR4
CXCL12
Tumor Invasion
Metastasis
Cytotoxicity
Imaging
Microscopy
Cellular Analysis
Phenotypic Assay
Label-free

BioTek Instruments, Inc.
P.O. Box 998, Highland Park,
Winooski, Vermont 05404-0998 USA
Phone: 888-451-5171
Outside the USA: 802-655-4740
Email: customercare@biotek.com
www.biotek.com
Copyright © 2015

treated microplates (Catalog No. 3674), used for performance of all ligand binding experiments, were donated by Corning Life Sciences (Kennebunk, ME).

Synergy™ Neo Multi-Mode Reader

Synergy Neo combines a filter- and monochromator-based detection system in one compact unit. The HTRF® certified reader uses a high performance xenon flash lamp and dual PMTs in the filter-based optics to simultaneously detect the assay's 665 nm and 620 nm fluorescent emissions when the excitation was set to 340 nm. The dual reagent injection capabilities enables kinetic analyses of known ligand and competitive binding. Tag-lite cells, in a volume of 10 μ L, were injected into the wells of the 384-well plates, already containing ligands and buffer, to begin each kinetic binding experiment.

Cytation™ 5 Cell Imaging Multi-Mode Reader

Cytation 5 is a modular multi-mode microplate reader that combines automated digital microscopy and microplate detection. Cytation 5 includes filter- and monochromator-based microplate reading; the microscopy module provides high resolution microscopy in fluorescence, brightfield, color brightfield and phase contrast. With special emphasis on live-cell assays, Cytation 5 features temperature control to 65 °C, CO₂/O₂ gas control and dual injectors for kinetic assays. Shaking and Gen5 software are also standard. The instrument was used to image spheroids, as well as individual cell invasion through the Matrigel matrix.

Tag-lite® Receptor Ligand Binding Assay




Figure 1. Tag-lite® Receptor Ligand Binding Assay Procedure.

The Tag-lite CXCR4 assay relies on a fully functional SNAP-tag fused CXCR4 receptor and fluorescently labeled ligand SDF1- α . Being homogeneous, the binding assay allows for binding events to be precisely recorded in time. The assay can be used to derive the kinetic binding parameters of unlabeled compounds by application of the Motulsky and Mahan equations.

Corning® Spheroid Microplates

Corning 96-well black, clear-bottom spheroid microplates (Catalog No. 4520) are coated with the Ultra Low Attachment surface which is a non-cytotoxic, and biologically inert covalently bonded hydrogel that prevents cell attachment. Novel well geometry aids spheroid formation in the center of each well. Each microplate contains an optically clear round bottom, which is ideal for cellular imaging, as well as a black, opaque body which prevents cross talk. The plates were used for performance of all 3D tumor invasion experiments.

Methods

Kinetic Binding

SDF1- α -d2 Association Rate Constant

The association rate constant of SDF1- α -d2 ($k_{on, SDF1-\alpha-d2}$) was determined by placing the CXCR4 expressing cells (5000 cells/well) in the presence of SDF1- α -d2 and measuring specific binding at various times thereafter. A total of seven [SDF1- α -d2] were tested, thus also allowing for calculation of the half saturation binding constant at equilibrium ($K_{d, SDF1-\alpha-d2}$).

SDF1- α -d2 Dissociation Rate Constant

The dissociation rate constant of SDF1- α -d2 ($k_{off, SDF1-\alpha-d2}$) was determined by adding to the previous mix a 100 fold molar excess of AMD3100, effectively blocking further binding of SDF1- α -d2 to its receptor.

Kinetic Characterization of Inhibitor Panel

The association and dissociation rate constant of the different CXCR4 inhibitors were tested following the kinetic competitive binding protocols described elsewhere. All experimental data were analyzed by using GraphPad Prism 5.0 (GraphPad Software, Inc., San Diego, CA).

Inhibition of Phenotypic Tumor Invasion

Cell Preparation and Spheroid Formation

MDA-MB-231 and fibroblast cells were harvested and diluted to a concentration of 5.0×10^4 cells/mL in complete medium. The two volumes were combined to create final concentrations of 2.5×10^4 cells/mL for each cell type. 100 μ L of cell suspension was then pipetted to appropriate wells. Following dispensing, the plate was placed at 37 °C/5% CO₂.

Image-based Spheroid Formation Monitoring

Spheroid formation was monitored every 24 hours. The plate was placed into the Cytation 5, previously set to 37 °C/5% CO₂ using Gen5 as well as a gas control module. Focusing was performed using the brightfield channel. The typical cell aggregation period was 48 hours.

Invasion Matrix Preparation

Upon spheroid formation completion, 70 µL of complete medium was removed from each well, washed with an equal volume of invasion medium (serum-free, phenol red-free medium), and the spheroid plate placed on ice in a refrigerator for 5 minutes to cool the wells. Corning® Matrigel® Matrix, Phenol Red-Free was then thawed on ice, while unlabeled SDF1- α was diluted to a concentration of 20 ng/mL in invasion medium. SDF1- α was further diluted to a concentration of 10 ng/mL in either invasion medium or Matrigel matrix. CXCR4 inhibitors were then titrated to 1x concentrations in both diluted invasion medium containing 10 ng/mL SDF1- α , and matrix containing 10 ng/mL SDF1- α .

With the plate still on ice, 70 µL of invasion medium with SDF1- α and inhibitor was added to each well. A 100 µL volume of matrix, containing SDF1- α and inhibitor, was then added as an overlay to each well. The plate was centrifuged at 300 x g for 5 minutes in a swinging bucket centrifuge that had been previously set to 4 °C for spheroid positioning, and then transferred to a 37 °C/5% CO₂ incubator for one hour to initiate gel formation.

Tumor Invasion Assay Performance

Using a 4x objective, exposure settings were optimized for the brightfield, as well as GFP and RFP fluorescent imaging channels. Following the optimization process, automated day 0 imaging was performed, and continued every 24 hours pursuant, to track tumor invasion. Cellular analysis was performed with the captured 4x images to track invasion in the presence or absence of CXCR4 inhibitors.

Results and Discussion

Drug-Target Residence Time Determination

Association Kinetics of SDF1- α -d2 Labeled Ligand

The final Drug-Target Residence Time value takes into account the observed on and off rates of the unlabeled inhibitors as well as the labeled SDF1- α -d2 ligand, and is computed by incorporation of the Motulsky and Mahan equation⁹. The first step to calculate the

final value was to perform an associative binding experiment using a concentration range of 0-100 nM of the d2 acceptor fluor labeled ligand. Binding was monitored kinetically over a period of 40 minutes.




Figure 2. Association binding graph of SDF1- α -d2. Observed associative binding curves calculated from HTRF ratios of wells containing SDF1- α -d2 ligand concentrations ranging from 0-100 nM. Non-specific binding values subtracted from total ratios to determine observed specific binding.

Binding increases over time until it plateaus after several minutes (Figure 2). The plateau in an association experiment depends on the concentration of labeled SDF1- α used. Higher plateaus will be obtained with higher concentrations. Fitting of the curves with Graph Pad Prism yields the observed association rate values for all concentrations tested or k_{obs} .

The K_d value of the labeled ligand was also determined by plotting the HTRF ratios generated after a binding equilibrium was reached with the different concentrations of ligand tested.




Figure 3. SDF1- α -d2 saturation binding curve. HTRF ratios generated upon the achievement of binding equilibrium of tested [SDF1- α -d2].

In a saturation binding experiment, increasing concentrations of labeled SDF1- α result in increased binding. Saturation is obtained when no further binding can be recorded. The ligand concentration that binds to half the receptor sites at equilibrium or K_d was 29 nM.

An assessment of whether the labeled SDF1- α ligand follows the Law of Mass action can also be carried out. If the system does follow the Law of Mass action then k_{obs} increases linearly with increasing concentrations of SDF1- α .




Figure 4. Law of Mass Action Determination. Plot of k_{obs} values calculated in relation to [SDF1- α -d2].

Due to the linear shape of the curve, and an R^2 value >0.9 , Law of Mass Action was proven for the labeled SDF1- α ligand. This allowed for the use of Graph Pad Prism software to derive association and dissociation rate constants from the linear regression line.

The rate constant values experimentally found or mathematically derived are summarized in Table 1. $k_{on, SDF1-\alpha-d2}$ and $k_{off, SDF1-\alpha-d2}$ were $0.001 \text{ nM}^{-1} \cdot \text{s}^{-1}$ and 0.04 s^{-1} , respectively.

| SDF1- α -d2 Kinetic Binding Characterization | | | | | | | |
|--------------------------------------------------------------------------------------------------|---------|---------|---------|---------|---------|---------|---------|
| [SDF1- α -d2] (nM) | 1.562 | 3.125 | 6.25 | 12.5 | 25 | 50 | 100 |
| $k_{off} (\text{s}^{-1})$ Experimental | 0.045 | 0.031 | 0.034 | 0.041 | 0.042 | 0.043 | 0.034 |
| $k_{obs} (\text{nM}^{-1} \cdot \text{s}^{-1})$ Experimental | 0.003 | 0.006 | 0.007 | 0.011 | 0.012 | 0.015 | 0.019 |
| $k_{on} (\text{nM}^{-1} \cdot \text{s}^{-1})$ calculated from $k_{obs} = k_{on} * [L] + k_{off}$ | 0.03033 | 0.01201 | 0.00647 | 0.00415 | 0.00217 | 0.00115 | 0.00054 |
| $k_{on} (\text{nM}^{-1} \cdot \text{s}^{-1})$ calculated from $K_d = k_{off} / k_{on}$ | 0.0010 | | | | | | |
| K_d (nM) Experimental | 29.3 | | | | | | |

Table 1. SDF1- α Kinetic Binding Characterization.

Association Kinetics of SDF1- α -d2 Labeled Ligand

In the theory developed by Motulsky and Mahan, an unlabeled competitor is co-incubated with a labeled ligand during a kinetic association experiment. Here, a single concentration of the SDF1- α -d2 ligand, 25 nM, was co-incubated with multiple concentrations of the unlabeled SDF1- α competitors in the presence of the CXCR4 expressing cells. Kinetic binding of the labeled ligand was then monitored over time.




Figure 5. Kinetics of Competitive Binding. Plot of specific binding HTRF ratios over time for the SDF1- α -d2 ligand when in the presence of 100, 10, or 1 nM concentrations of (A.) AMD 3100, (B.) AMD 3465, or (C.) IT1t.

From the curve fitting of the observed SDF1- α -d2 kinetic binding, and incorporation of the Law of Mass Action linear regression line, k_{off} (Min^{-1}) values were then calculated. Final residence time (R) values could then be determined using the following formula:

$$R = 1/k_{\text{off}}$$

Therefore, molecules having a lower k_{off} rate reside at the target receptor for longer periods of time.

| | k_{off} (Min^{-1}) | R (Min) | K_i (M) |
|------------------------------------|----------------------------------------|---------|-----------|
| SDF1-α-d2 | 2.28 | 0.44 | 2.40E-08 |
| AMD3465 | 0.003044 | 328.52 | 1.44E-10 |
| IT1t | 0.02844 | 35.16 | 1.93E-10 |
| AMD3100 | 5.496 | 0.18 | 1.56E-08 |
| CTCE9908 | 8.692 | 0.12 | 3.51E-06 |
| FC131 | ND | ND | 4.52E-09 |

Table 2. SDF1- α Competitor Dissociation Rate and Residence Time Values.

From the shape of the curves in Figure 5, and a comparison of the residence time values generated for the labeled ligand and unlabeled competitors (Table 2), qualitative and quantitative assumptions regarding the various competitors can then be made. First, if the competitor dissociates faster from its target than the ligand (smaller R value), such as is seen with AMD 3100 (Figure 5A), the specific binding of the ligand will slowly and monotonically approach its equilibrium in time. However, when the competitor dissociates slower (larger R value), the association curve of the ligand consists of two phases, starting with a typical "overshoot" and then a decline until a new equilibrium is reached. Competitors whose residence times are greater than that of the SDF1- α -d2 ligand, such as AMD 3465 and IT1t (Figure 5B and C), may then exhibit a stronger inhibitory response when used in the confirmatory phenotypic 3D tumor invasion assay.

Interruption of Invasion via SDF1- α Ligand Binding Inhibition

As stated previously, interruption of the interaction between CXCR4 and its known ligand, SDF1- α , impairs metastasis of breast cancer and cell migration². Therefore, a phenotypic assessment of the CXCR4 inhibitor panel was then performed to determine whether changes in the level of tumor migration could be detected, and more importantly, if compounds exhibiting longer residence times compared to SDF1- α -d2 exhibited a higher inhibitory effect on migration through the 3D matrix. MDA-MB-231 breast adenocarcinoma cells, co-cultured with human dermal fibroblasts, were used as the *in vitro* tumor model. This breast cancer cell line has been previously shown to express the CXCR4 receptor¹⁰.




Figure 6. Image-based Monitoring of MDA-MB-231/Fibroblast Tumor Invasion. Overlaid brightfield and fluorescent images captured using a 4x objective, after a 0 and 5 day incubation period with AMD 3465, IT1t, and CTCE 9908. Imaging channel representation: Brightfield – Total cells and invadopodia; GFP – MDA-MB-231 cells; RFP – Fibroblasts.




Figure 7. Quantification of Invasive Tumor Area. 4x overlaid images captured following 5 day (A.) 100 and (B.) 0 μ M IT1t incubation with tumoroids. Object masks automatically drawn by Gen5 using the following criteria: Threshold: 5000 RFU; Min. Object Size: 400 μ m; Max. Object Size: 1500 μ m; Image Smoothing Strength: 0; Background Flattening Size: Auto.

Cellular analysis is performed with the Cytation 5 using the brightfield signal to quantify the extent of invasion. Minimum and maximum object sizes, as well as brightfield threshold values are set such that a precise object mask is automatically drawn around each tumoroid in its entirety (Figure 7A and B). The same criteria are used for all images evaluated during the experiment. This allows for a quantitative comparison of the area covered within each object mask to be completed.




Figure 8. Tumor Invasion Inhibition Determination. Graphs of individual tumoroid areas on day 0, and subsequent to five day invasion period in the presence of inhibitor concentrations.

The 4x images displayed (Figure 6), as well as the graphs in Figure 8, demonstrating total tumoroid area coverage before and after the incubation period illustrate the ability of CXCR4 inhibitors to interrupt tumor invasion consistent with the previously determined residence time. AMD 3465 and IT1t, which exhibit a residence time longer than SDF1- α -d2, effectively minimize tumor invasion in a dose dependent manner. The decrease in MDA-MB-231 GFP and fibroblast RFP expression exhibited after a 5 day 100 μ M IT1t incubation, also seen after a 7 day AMD 3465 incubation of the same concentration (data not shown), may also indicate the chronic cytotoxic effects that elevated dosing of these compounds can have on both cancer and stromal cells. All other compounds show little to no effect on the ability of the tumoroid to migrate through the 3D matrix. While AMD 3465 and IT1t display the same sub-nanomolar potency, AMD3465 prevails as a CXCR4 inhibitor due to its greater residence time.

Conclusions

The Tag-lite CXCR4 ligand binding assay provides a simple, yet robust cell-based approach to determine kinetic binding of known receptor ligands, as well as competitive binding of test molecules. The simultaneous dual emission capture and injection capabilities of the Synergy Neo allow accurate calculations of kinetic association and dissociation rates to be made when used in conjunction with the Tag-lite® assay. Corning Spheroid Microplates then provide an easy-to-use, consistent method to perform spheroid aggregation and confirmatory 3D tumor invasion assays. Imaging of spheroid formation, as well as invading structures can be performed by the Cytation™ 5 using brightfield or fluorescent channels to easily track tumoroid invasion. The flexible cellular analysis capacity of the Gen5™ Data Analysis Software also allows for accurate assessment of 3D tumor invasion during the entire incubation period. The combination of assay chemistry, cell model, kinetic microplate and image-based monitoring, in addition to cellular analysis provide an ideal method to better understand the target-based and phenotypic effects of potential inhibitors of tumor invasion and metastasis.

References

1. Saxe, Charles. 'Unlocking The Mysteries Of Metastasis'. *ExpertVoices* 2013. <http://www.cancer.org/cancer/news/expertvoices/post/2013/01/23/unlocking-the-mysteries-of-metastasis.aspx>. Accessed 16 Mar. 2015.
2. Müller, A., Homey, B., Soto, H., Ge, N., Catron, D., Buchanan, M., McClanahan, T., Mruphy, E., Yuan, W., Wagner, S., Barrera, J., Mohar, A., Verástegui, E., Zlotnik, A. Involvement of chemokine receptors in breast cancer metastasis. *Nature*. **2001**, 410, 50-56.
3. Swinney, D. Biochemical mechanisms of drug action: what does it take for success? *Nat Rev Drug Discov*. **2004**, 3, 801-808.
4. Copeland, R., Pompliano, D., Meek, T. Drug-target residence time and its implications for lead optimization. *Nat Rev Drug Discov*. **2006**, 5, 730-739.
5. Tummino, P., Copeland, R. Residence time of receptor-ligand complexes and its effect on biological function. *Biochemistry*. **2008**, 47, 5481-5492.
6. Zhang, R., Monsma, F. The importance of drug-target residence time. *Curr Opin Drug Discov Devel*. **2009**, 12, 488-496.
7. Mao, Y., Keller, E., Garfield, D., Shen, K., Wang, J. Stromal cells in tumor microenvironment and breast cancer. *Cancer Metast Rev*. **2013**, 32, 303-315.
8. Kamath, L., Meydani, A., Foss, F., Kuliopulos, A. Signaling from protease-activated receptor-1 inhibits migration and invasion of breast cancer cells. *Cancer Res*. **2001**, 61, 5933-5940.
9. Motulsky, H., Mahan, L. The kinetics of competitive radioligand binding predicted by the law of mass action. *Mol Pharmacol*. **1984**, 25, 1-9.
10. Sun, Y., Mao, X, Fan, C, Liu, C., Guo, A., Guan, S., Jin, Q., Li, B., Yao, F., Jin, F. CXCL12-CXCR4 axis promotes the natural selection of breast cancer cell metastasis. *Tumor Biol*. **2014**, 35, 7765-7773.