


## Facilitated Process Scalability Using a Novel Raman Spectroscopy Integration

Michael Sibley<sup>1</sup>, Angus Woodhams<sup>1</sup>, Marek Hoehse<sup>2</sup>, Svea Cheeseman<sup>1\*</sup>  
<sup>1</sup>Sartorius Royston, UK. <sup>2</sup>Sartorius Göttingen, Germany. \*svea.cheeseman@sartorius.com

### Aims

- Demonstrate hardware comparability of the BioPAT® Spectro interfaces in Ambr® and Flexsafe STR® and existing immersion probes
- Demonstrate model scalability of a cell culture process from Ambr® to the production size Flexsafe STR® bags

### Key Features of BioPAT® Spectro in Ambr® and Flexsafe STR®



### BioPAT® Spectro Key Features


- Scalable Raman Spectroscopy platform from 15 mL to 2000L
- Ambr® Analysis Module with replaceable Raman flow cell
- Compatibility: Ambr® 15 Cell Culture, Ambr® 250 High Throughput
- Fully automated data acquisition, alignment and model application in Ambr®, with automated spiking
- Fully integrated and sterilized single-use port in Flexsafe STR® 50 – 2000 L
- Port qualified for GMP manufacturing, shields measurement from ambient light and decouples the measurement from process conditions such as gassing and stirring
- Raman spectroscopy: Compatibility with Raman spectrometers from Kaiser Optical Systems and Tornado Spectral Systems
- Facilitated model transfer from Ambr® to Biostat STR® and vice versa due to identical optical design and the use of the same fiber optic probe



### Hardware Scalability Tested in a Lab Set-Up

#### Model building using the Ambr® flow cell

- Ambr® flow cell with four different optical windows
- Raman spectra were acquired for each combination with glucose concentrations between 0.5 and 10 g/L
- Four corresponding glucose models were built using SIMCA®




#### Model transfer from the Ambr® flow cell to the Flexsafe STR® port and an immersion probe

- Single-use ports of the Flexsafe STR® with the same four optical windows from above
- Raman spectra were acquired for each combination with glucose concentrations between 0.5 and 10 g/L
- A conventional immersion probe used to acquire spectra in a beaker with 0.5-10g/L glucose
- Glucose multi-variate models generated with the Ambr® flow cells were used to predict the glucose concentrations



#### Model transfer results


- Ambr® models led to a consistent and precise prediction of the glucose concentrations in the Flexsafe STR® ports and the immersion probe
- The average prediction error across all predictions was 0.23 g/L and is comparable to the literature.



### Process Scalability Tested in a Cell Culture Process

#### Model building using the Ambr® 250 High Throughput

- N = 16 Ambr® 250 High Throughput bioreactors
- Sartorius' Cellca2 Process CHO producing mAb
- Integrated Raman and Nova FLEX2 analyzers
- 200 data points assayed by Flex2 then Raman, of which ~35 were spiked with 7.1 or 33.3 g/L glucose
- Data was automatically collated in the Ambr® software and manually exported to SIMCA®
- SIMCA 16 software was used for multivariate model building




#### Model transfer from Ambr® 250 to 200L Biostat STR®

- Cellca 2 process was run in a 200L Flexsafe STR® bag, at a different Sartorius site using a spectrometer and fiber optic probe with different serial numbers
- BioPAT® Spectro SU port was used, which comes integrated and sterilized inside the bag
- The Ambr® multi-variate model was applied to the spectra without further modifications to predict glucose

#### Model transfer results

- The Ambr® model on the Biostat STR® data predicted the glucose profile and identifies feeding events
- Predictions agreed with offline measured glucose concentration
- Improved results are to be expected when updating the Ambr® Raman model with spectra acquired in at least one Biostat STR® run, which we recommend



### Conclusions

- Raman models transfer between the hardware used to acquire Raman spectra in the Ambr® and Biostat STR® systems in a laboratory set-up
- The BioPAT® Spectro optical interfaces are comparable to conventional immersion probes and models can be transferred between them in a laboratory set-up
- The Raman glucose model for the Cellca2 process in the Ambr® 250 High Throughput acquired by the BioPAT® Spectro integration is readily transferable to the same process in the Flexsafe STR® bag using the BioPAT® Spectro SU port
- This indicates that models built for other processes and analytes are also likely to scale well in the BioPAT® Spectro platform