


# In Vitro Test Methodologies for Characterizing Bioavailability Enhancing Formulations

**Aaron Stewart**

***Associate Principal Scientist, Global R&D***

©2020 Lonza Pharma & Biotech | Webinar | In Vitro Test Methodologies for Characterizing Bioavailability Enhancing Formulations


# Talk Outline




Strategies to Address  
Low Bioavailability




In Vitro Performance Tools  
Methodologies  
How to select?



Case Study #1  
Itraconazole




Case Study #2  
Belinostat



Summary and  
Take Away Points

# The Ongoing Issue

## Low solubility continues to plague development pipelines



Reference:

*Drug-Like Property Concepts in Pharmaceutical Design, Di, Li; Kerns, Edward H.; Carter, Guy T. In Current Pharmaceutical Design, Volume 15, Number 19, 2009, pp. 2184-2194(11)*


# Putting Low Drug Solubility in Context

## Primary obstacles to low solubility

### Solid-state



“Brick dust” compounds




- Very strong crystalline lattice energy **(1)** is limiting to solubility
- Address by bypassing crystalline state

### Solvation



“Greaseball” compounds



- Very low affinity for water **(3)** is limiting to solubility
- Address by bypassing crystalline state and/or changing the nature of the solvent – “like dissolves like”

It is important to know the solubility performance attributes for your compound as it will aid selecting the right technology

# Many Enabling Technologies Are Available

## Solid-State Alteration: Form, Particle Size


- Polymorphs
- Amorphous Solid dispersions
  - SDD, HME
  - Lyophiles
  - Drug/polymer nanoparticles
  - Mesoporous carriers
- Nanocrystals (200 to 800 nm)
- Nanocrystals (<50 nm)

## New crystalline compound

- Cocrystals
- Salts

## Solvation, Complexation


- Cosolvents
- Surfactants
- Cyclodextrins
- Lipids:
  - Oils
  - SEDDS/SMEDDS
  - Solid lipid pellets



H.D. Williams et al. "Strategies to Address Low Solubility in Discovery and Development," Pharmacol Rev , 65(2013)315-499

# Which Technology To Select for Which Compound?

Why it is an important question....



## Matching your compound with the wrong technology

Challenging development

Increase development time and costs

High risk of failure



## Exploring multiple technologies in parallel



Complex development

Higher resource requirement (API, costs)

Does not guarantee success

# Compound Qualification and Technology Mapping

A science-based approach to selecting a technology for your compound



- Broad, in-depth review considers physical and biological barriers to drug absorption
- Drawbacks: Qualitative assessment. Does not factor in experience and technology precedence

- Leverages large compound *in vivo* datasets and based on data collected during preformulation
- Drawbacks: Data does not include head-to-head comparisons. Would another technology have been better?

# Representative Technology Options Relative to Problem Statement

Technology	Dissolution Rate Limited	Permeability Limited	Solubility and Permeability Limited
Salt, Polymorph, Cocrystal	X		X
Nanocrystals	XX		
Amorphous	XX		X
Lipids	XX	X	XX

# In Vitro Bioperformance Tools

Available tools...

How to select?



# *In vitro* & *in silico* tools can be used to assess the rate determining step(s) to absorption

## Basic property prediction



ChemAxon Ltd

- HBD, HBA, PSA
- Acid/base/neutral
- pK<sub>a</sub>
- LogP/LogD
- pH-solubility
- T<sub>m</sub>
- Provisional BCS

## Dimensionless numbers/simple calculations

- Dose Number
- Dissolution Number
- Permeation Number

## *In vitro* assessment

- pH-solubility
- LogP/LogD
- T<sub>m</sub> & T<sub>g</sub>
- Micelle partition coeff.
- Caco2 permeability
- PAMPA permeability

# Dimensionless Numbers

- Key Inputs:
- Crystalline Aqueous Solubility
- Amorphous Aqueous Solubility
- Projected Dose (or range)
- Animal physiology
- Log P/Micelle Partition
- Effective Permeability

Outputs:

Dose Number

$$Do = \frac{Dose/Vol}{C_s}$$



How many GI fluid volumes are required to dissolve the dose ?

Permeation Number

$$Pn = k_{abs} \cdot t_{abs}$$



How many times can the drug permeate over the course of GI transit ?

Dissolution Number

$$Dn = k_{diss} \cdot t_{abs}$$



How many times can the drug dissolve over the course of GI transit ?

# Fraction Absorbed Classification System (FaCS)

## Three Limiting Cases

Dose Number

$$D_o = \frac{\text{Dose}/\text{Vol}}{C_s}$$


Permeation Number

$$P_n = k_{abs} \cdot t_{abs}$$

Dissolution Number

$$D_n = k_{diss} \cdot t_{abs}$$

Case 1: Dissolution Rate Limited (DRL)



Cases where this occurs:

- High permeability relative to dose and solubility
- Dissolution rate is slow

# Fraction Absorbed Classification System (FaCS)

## Three Limiting Cases

Dose Number

$$Do = \frac{Dose / Vol}{Cs}$$


Permeation Number

$$Pn = k_{abs} \cdot t_{abs}$$


Dissolution Number

$$Dn = k_{diss} \cdot t_{abs}$$

Case 1: Dissolution Rate Limited (DRL)



Case 2: Permeability limited (PL)



Cases where this occurs:

- High permeability relative to dose and solubility
- Dissolution rate is slow

Cases where this occurs:

- Dissolution rate is fast relative to permeability
- Dose is low relative to solubility

# Fraction Absorbed Classification System (FaCS)

## Three Limiting Cases

Dose Number

$$D_o = \frac{\text{Dose}/\text{Vol}}{C_s}$$


Permeation Number

$$P_n = k_{abs} \cdot t_{abs}$$


Dissolution Number

$$D_n = k_{diss} \cdot t_{abs}$$


Case 1: Dissolution Rate Limited (DRL)



Case 2: Permeability limited (PL)



Case 3: Solubility-permeability limited (SL)



Cases where this occurs:

- High permeability relative to dose and solubility
- Dissolution rate is slow

Cases where this occurs:

- Dissolution rate is fast relative to permeability
- Dose is low relative to solubility

Cases where this occurs:

- Low permeability relative to dose and solubility
- Dose is high relative to solubility

# How do we select and design the appropriate *in vitro* test?

## 1. Predict *in vivo* problem statement

Problem Statement	Formulation mechanism
Precipitation	Supersaturating formulations
	pH driven supersaturation
Dissolution Rate Limited	Initial particle size
	Dissolution rate changes over time
Solubility-Permeability Limited	ABL limited
	Epithelial limited

## 2. Select *in vitro* dissolution apparatus

<i>In Vitro</i> toolkit	
Non-sink Single medium	Dissolution – permeation
Non-sink pH shift	Controlled Transfer
Sink single medium	Dissolution - permeation
Sink pH shift	Controlled Transfer
Non-sink Single medium	Dissolution - permeation
Non-sink pH shift	

## 3. Choose *in vitro* test parameters

<i>In Vitro</i> parameters	
Fluid composition, volume	A/V, volume(s)
	Fluid transfer rates
Fluid composition	A/V, volume(s)
	Fluid transfer rates
Fluid composition, volume	A/V, volume(s)

# Multiple Problem Statement-specific Bioperformance *In Vitro* Tools Using Fiber Optics

Amorphous Solubility



- Amorphous “solubility”
- Precipitation risk
- Polymer selection
- Drug/polymer interaction

Dissolution




- Dissolution rate
- Precipitation rate
- Maximum apparent concentration
- Speciation

Flux



- Clean measurement of “effective” concentration
- Able to properly account for micelle, colloid, and particle contribution to boundary layer diffusion and dissolution rate

Controlled Transfer





- Dissolution rate
- Precipitation rate vs. emptying rate
- Gastric precipitation
- “Book-end” for formulation performance

# Amorphous Solubility and Polymer Screening

## Key Outputs:

- Maximum solubility attainable from amorphous formulation
- Precipitation risk of compound
- Identify lead sustaining polymer(s) for SDD
- Potential to identify strong drug/polymer interaction



Amenable to very low API quantities of 20 mg or less using a small volume setup



Murdande *et al.* *J Pharm Sci.* 99 (2010), 1254  
Murdande *et al.* *Pharm Res.* 27 (2010), 2704  
Almeida e Sousa *et al.* *Mol. Pharmaceutics.* 12 (2015), 484  
Ilievare and Taylor. *Cryst. Growth Des.* 13 (2013), 1497

# Pion UV Probe Dissolution

## Key Outputs:

- Dissolution rate/extent
- Precipitation (gastric or intestinal)
- Speciation customized for specific API (e.g unbound free drug, micelle bound drug, colloids)


*Example gastric to intestinal pH transfer test*



# Membrane Flux

## Key Outputs:


- Measures “effective concentration” for species essential for membrane permeation
- Can account for species contribution to boundary layer diffusion and dissolution rate




Concentration gradient example above demonstrates an Aq. boundary layer-limited formulation

### Resistances in series

- Dissolution limited
- Aq. Boundary layer limited
- Membrane limited



# Measured Flux is a Sum of Resistances



## Dissolution Rate

- Dose
- Volume
- Solubility
- Particle size

## Boundary Layer Diffusion


- Unbound drug
- Micelle-bound drug
- Colloidal drug

## Membrane Diffusion

- Unbound neutral drug

# Membrane Flux

## How Data Inform Formulation Decisions




Route highlighted in green typical example of a low aqueous solubility, highly lipophilic BCS Class II compound


# Controlled Transfer Dissolution

## Key Outputs:

- Biorelevant dissolution rate
- Precipitation rate at biorelevant mass/volume transfer
- Gastric precipitation
- Dosage form disintegration



## Solid Dosage Form Scale



## Intermediate Scale





# Digestion test for lipid based formulations


Digestion tests are used to identify those formulations at low risk of showing drug precipitation in the intestine

Provides a performance-based means to select lead formulations

## Measure digestion by titration



## Sample and measure drug fate by HPLC




Recent success has been demonstrated incorporating an absorption compartment into the digestion test: Alvebratt, C., Keemink, J., Edueng, K., & Cheung, O. (2020). *European Journal of Pharmaceutics and Biopharmaceutics*. <https://doi.org/10.1016/j.ejpb.2020.01.010>


# Case Study 1: Itraconazole



# Case study - Amorphous spray dried dispersions (SDDs) of Itraconazole (ITZ) dosed to rats



**Hydroxypropyl Methylcellulose Acetate Succinate (HPMCAS)**




**25% active  
Hydrophilic SDD  
Affinisol 716HP**



or

**25% active  
Hydrophobic SDD  
Affinisol 126HP**



## Formulations dosed to rats

Sprague-Dawley (n=6), fasted

Dose: 50 mg/kg

Dosing vehicle: 0.5% Methocel A4M

in H<sub>2</sub>O

Dosing route: oral gavage

## Key ITZ attributes:

- Itraconazole has exceedingly low aqueous solubility even in the amorphous state (ca. 0.1 µg/mL).
- High lipophilicity and neutral charge state at intestinal pH drives very low solubility but high lipid membrane permeability, resulting in aqueous boundary layer limited flux *in vitro*.

# Dimensionless numbers can predict impact of solubility, permeability or dissolution rate in vivo for itraconazole

**BCS** Ref: Amidon, G.L., et al. *Pharm Res.* (1995), 12 (3), 413-420

**FaCS** Ref: Sugano, K., et al., *J Pharm Sci.* (2015), 104, 2777-2788

## Dose Number


$$Do = \frac{Dose / Vol}{C_s}$$

## Dissolution Number


$$Dn = k_{diss} \cdot t_{abs}$$

## Permeation Number


$$Pn = k_{abs} \cdot t_{abs}$$




## Dissolution-limited





## Permeability-limited



# Itraconazole is highly solubilized in micelles and colloids





# Material sparing in vitro membrane flux test can assess solubility-permeability limited absorption



Itraconazole = ABL Limited Flux

# Hydrophilic SDD has the highest flux *in vitro*



No.	Formulation	Dispersion polymer	Flux (µg/min/cm <sup>2</sup> )	Colloid (µg/ml)
■	25% ITZ/75% HPMCAS SDD	AFFINISOL 716HP	1.18	602
●	25% ITZ/75% HPMCAS SDD	AFFINISOL 126HP	0.85	150
◆	Sporanox <sup>®</sup> spray layered dispersion	HPMC	0.53	0

All formulations have the same unbound (0.1 µg/mL) and micelle-bound (20 µg/mL) ITZ concentrations and only differ in the concentration of colloidal drug species. Difference in flux is driven by the nano-sized colloidal species.

# ABL limited diffusion in the membrane flux assay can be described by a steady state diffusion model


$$j = \frac{D_{eff}}{h_{ABL}} (c_{u,m,c})$$

$$D_{eff} = D_u \cdot f_u + D_m \cdot f_m + D_c \cdot f_c$$

$C_{u,m,c}$  is the sum of all species:



$D_c = 4 \times 10^{-8} \text{ cm}^2/\text{s}$   
(200 nm)




$D_m = 1 \times 10^{-6} \text{ cm}^2/\text{s}$   
(7 nm)



$D_u = 5 \times 10^{-6} \text{ cm}^2/\text{s}$

## Key assumptions



- Pseudo-steady state
- $D_u$ ,  $D_m$  &  $D_c$  concentration independent
- Drug that is unbound, micelle-bound, or in colloids contribute to  $D_{eff}$  based on size and abundance
- Minimal transport due to convection
- Well mixed solutions
- Constant ABL thickness



Model supports *in vitro* measurements made in three different media: blank PBS buffer, 6.7 mM SIF, 27 mM SIF.

Dose for all flux measurements was 1000  $\mu\text{g/mL}$  ITZ

# Hydrophilic SDD shows the fastest absorption in rats – rank orders with *in vitro* performance



# Itraconazole Case Study


## Conclusions

- Identified unique drug speciation from ITZ:Affinisol SDDs compared to commercial formulation Sporanox
- Evaluated contributions of these species to in vitro flux based on ABL limited diffusion
- Described contributions of drug species mathematically
- Demonstrated the impact *in vivo*, showing absorption rate trends with *in vitro* flux.
  
- Key in vitro performance tool: **membrane flux**

# Case Study 2: Belinostat



# Case study - SDDs of belinostat dosed to dogs



**Belinostat**

BCS II/IV


$pK_a = \geq 8$  (acidic)

$\text{LogP} < 2$

HPMCAS (weakly acidic)



25% active  
HPMCAS-M SDD



+ Polyvinylpyrrolidone (neutral)




25% active  
PVP K30 SDD



Polyvinylpyrrolidone Vinyl Acetate (neutral)



25% active  
PVP VA64 SDD



## SDDs dosed to beagle dogs

(n=4), fasted

Dose: 50 mg

Dosing vehicle: 0.5% Methocel

A4M in H<sub>2</sub>O, 15 ml water rinse

Key belinostat attributes:

- High amorphous solubility in biorelevant media (>500 µg/mL).
- Amorphous solubility is impacted by the presence of polymer.
- Dissolution rate is a key driver for absorption and differs depending on SDD formulation and testing method.



# Belinostat apparent amorphous solubility depends upon dispersion polymer type

## Belinostat

BCS II/IV

$pK_a = \geq 8$  (acidic)

$\text{LogP} < 2$



Amorphous solubility is defined as the onset of amorphous liquid-liquid phase separation. Presence of polymer influences the LLPS concentration.

# Dimensionless numbers can predict impact of solubility, permeability or dissolution rate *in vivo* for belinostat


**BCS** Ref: Amidon, G.L., et al. *Pharm Res.* (1995), 12 (3), 413-420

**FaCS** Ref: Sugano, K., et al., *J Pharm Sci.* (2015), 104, 2777-2788

## Dose Number

$$Do = \frac{Dose / Vol}{C_s}$$


## Solubility-permeability limited



## Dissolution Number

$$Dn = k_{diss} \cdot t_{abs}$$


## Dissolution-limited



## Permeation Number

$$Pn = k_{abs} \cdot t_{abs}$$

## Permeability-limited



# Evaluate belinostat dissolution performance using pH transfer test versus single medium test



*In situ* fiber optic detection

**Gastric transfer test**  
(pH 2 SGF → 6.5, 6.7 mM SIF)

Add Concentrated SIF solution at t = 30 min

**pH 2 SGF**      **pH 6.5**  
**6.7 mM SIF**

Non-sink Dose: 1000 µg/mL in SGF

	<i>In vitro</i> Gastric	<i>In vitro</i> Intestinal
HPMCAS-M SDD	1.3	0.4
PVP K30 SDD	1.4	0.4
PVP VA64 SDD	3.3	1.0

**Intestinal pH test**  
(pH 6.5, 6.7 mM SIF)

**pH 6.5**  
**6.7 mM SIF**

Non-sink Dose: 2000 µg/mL in SIF

	<i>In vitro</i> Intestinal
HPMCAS-M SDD	1.5
PVP K30 SDD	1.7
PVP VA64 SDD	4.0

***In vivo***

source: daviddarling.info


Assumes:

- Fasted state
- 50 mL gastric volume
- 50 mL intestinal volume


	<i>In vivo</i> Gastric	<i>In vivo</i> Intestinal
HPMCAS-M SDD	1.3	0.8
PVP K30 SDD	1.4	0.8
PVP VA64 SDD	3.3	2.0

Test design should be optimized towards the anticipated dose number and conditions *in vivo*.

# Relative extents of dissolution between SDDs depends upon dissolution medium composition





# Gastric → intestinal transfer test better rank orders SDDs with respect to *in vivo* performance in dogs






Sequential exposure to SGF and SIF at a more relevant dose/volume/solubility (dose number) is a better indicator for rank-ordering *in vivo* exposure from each SDD.

# Using amorphous solubility and dissolution data as key inputs to absorption model supports hypothesis of dissolution rate limited absorption

*In vitro*  
inputs to  
model



*In silico*  
predictions



# Belinostat Case Study

## Conclusions

- Amorphous solubility of belinostat depends on polymer type.
- SGF/SIF transfer test a better indicator of *in vivo* performance.
- Used *in vitro* inputs to describe blood plasma profiles via absorption modeling.
- Rate-determining step to absorption: dissolution rate and extent achieved in the stomach prior to transit down the GI tract.
  
- Key *in vitro* performance tool: **solvent shift amorphous solubility + fiber optic probe dissolution**

- The need for enabling technologies for improving oral bioavailability is not going away, and is likely going to increase moving forward.
- There is an industry driver to continue to develop an understanding of how to characterize and use in vitro data for supersaturating BAE formulations.
- With careful considerations for in vitro test methodology and design, we can bridge the gap from every compound being a “research project” to a platform-minded in vitro approach for BAE formulations in general.
- This is enabled through mechanistic understanding of bioperformance from BAE formulations independent of platform technology.
  - What is the rate limiting step to absorption?
  - What formulations characteristics will address this rate limiting step?
  - What in vitro tools do we have that can evaluate BAE formulation based characteristics?
  - In vivo relevance?

# Acknowledgements

## Lonza Pharma & Biotech colleagues

- Michael Grass, PhD
- Deanna Mudie, PhD
- Ian Yates
- Michael Morgen, PhD
- David Vodak, PhD
- David Lyon, PhD

The information contained herein is believed to be correct. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information. Lonza disclaims any liability for the use of this presentation and the use of the information contained herein is at your own risk.

All trademarks belong to Lonza or its affiliates or to their respective third-party owners.

Lonza 2020 ©

[Contact Us](#)